

Basi² Improved adjustability & Rigidity combined with a modern look!



Provides an even broader range of clinical solutions Breezy





Basi

Key Product Features:

Minimize Inventory:

• 1 Chair - Multiple Depths

Adjust Seat Depth by sliding cross brace 16" to 18" deep

• 1 Chair - Enhanced wheel adjustability to suit more clients

2

- 4 vertical positions, 2 active and 2 passive horizontal wheel positions
- Self Propel or Transit wheels available
- Tension adjustable backrest now extendable from 41 to 46 cm

Upgraded Features:









Upholstery

New Options:



arm rests



Adjustable Handles

Axle Plate

Caster Housing



Wheel Lock **Technical Specifications:**

Max. User Weight:	125 kgs
Chair weight	14.9 kgs
Overall width	SW + 19cm / SW +20cm (drumbrake version)
Seat width	14, 16, 18 & 20"
Seat Depth	16 - 18" - adjustable
Seat Height	48 cm (adj. 45.5 - 51 cm)
Back Height	41 - 46 cm
Armrest	20 - 30 cm (covers short & full length pad in one)
Armrest height adjustable	20.5 - 23 cm
Height adjustable push handles	+12 cm
Hangers:	70° (lower leg length adjustable from 35 - 50 cm) Angle adjustable footplate (+/- 15°)
Colour:	Robust hammerite - makes scratches less visible Red - on selected SKU's
Crash Test Approved	Yes
TGA#	145834
Warranty	Limited Lifetime



© Sunrise Medical. Breezy Basix 2. 2015. Rev B. AU

Sunrise Medical Australia Pty Ltd 6 Healey Circuit, Huntingwood, NSW 2148 Phone: (02) 9678 6600 Fax: (02) 9678 6655 Email: enquiries@sunrisemedical.com.au www.sunrisemedical.com.au