

Minos Plena

The budget alternative in aluminium

Minos Plena is an aluminium frame wheelchair with probably the highest quality in its segment.



For the most demanding users Minos Plena offers an adjustable fork and axle plate enabling the seat angle to be adjusted to fit the needs of the user.

Crash Tested

Minos Plena has been tested and approved for use in motor vehicles according to

ISO 7176
19:2001



MINOS

Think Global

New frame

Now you can switch from self-propelled to attendant-propelled in a few minutes with the new frame.



New adjustable in height armrest

You do not need any tool to quickly adjust the height of the armrest pad.

(Optional feature)



Dimensions

	Minos Plena	Minos Plena Transit
Seat width	36, 40, 44, 48, 51 cm	36, 40, 44, 48, 51 cm
Seat height	45-50 cm adjust.	45-50 cm adjust.
Seat depth	42 cm	42 cm
Backrest height	47 cm	47 cm
Armrest height	20-27 cm adjust.	20-27 cm adjust.
Total width	59, 63, 67, 71, 74 cm	50, 54, 58, 62, 65 cm
Total height	93 cm	95 cm
Total length	102 cm - 72 cm*	98 cm - 75 cm*
Folded width	31 cm	24 cm
Weight	14,8 kg *-**	13,1 kg *-**
Max. user's weight	125 kg ***	125 kg ***

* Without legrests.

** Weight without legrests. It can vary depending on accessories and options.

*** Minos Plena 51 (with double crossbar) Max. user's weight 160 kg.

Minos Global Accesories Compatible


The accessories and options of the Minos Global wheelchair are compatible with the Minos Plena.



Lifetime Guarantee

We are so convinced of the quality of the Minos Plena that we offer a lifetime guarantee for the crossbar and side frames.

Distributed by:



MINOS

Think Global

Pol. Ind. Sangróniz
Iberre Kalea, 7
48150 SONDIKA (España)

TEL. +34 94 453 32 30
FAX +34 94 453 30 01

E-MAIL pedidos@minos.es
WEB www.minos.es