## **J** AURORA







## FEATURES & BENEFITS

- Loop handlebars offer comfortable finger or thumb throttle operation
- 4 wheel adjustable suspension system allowing for more rider control and comfort
- Large 12" wheels which provide greater ground clearance
- Larger 4 pole motor and larger controller delivering much more power for climbing steep driveways or hills
- Separate emergency braking system for added security
- Built-in half speed safety switch which operates automatically when turning to provide extra safety

| Range           | 40-50 km                   |
|-----------------|----------------------------|
| Motor           | DC 24V 1750W 4 Pole        |
| Controller      | Rhino 160A                 |
| Battery         | 12V / 60-85 a/h x2 pieces  |
| Charger         | 24V/8amp off-board         |
| Gradient        | 10°                        |
| Front Wheels    | 12" pneumatic              |
| Rear Wheels     | 12" pneumatic              |
| Weight capacity | 200kg                      |
| Brake           | Intelligent, regenerative, |
|                 | electromagnetic brake      |
| Length          | 141cm (3 wheel)            |
|                 | 146cm (4 wheel)            |
| Width           | 65.5cm                     |