

Designed for showering over the bath, the lightweight base frame fits the majority of bath tubs.

The moulded plastic seat includes drain holes, comfortable contoured back rest and side support arms for added security.

The seat swivels 360° in both directions and includes a locking lever for set positioning and maximum user safety. The powder coated metal frame includes non-slip stoppers for maximum grip on the bath edge.

When using the Swivel Bath Seat, do not use any bath salts or other corrosive additive in the bath water. Doing so may cause corrosion which is not covered under the product warranty.





THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.

CODE	BL1522
CAPACITY	110 kg
OVERALL DEPTH	530 mm
OVERALL WIDTH	730 mm
OVERALL HEIGHT	500 mm
SEAT WIDTH	430 mm
SEAT DEPTH	430 mm
SEAT HEIGHT	355 mm
BETWEEN ARMRESTS	500 mm
CONFORMS TO	ARTG 176831 AS 4069-2013
WARRANTY	1year