

# Female Urinal Standard

Sturdy handle

Free of sharp edges,  
inside and out



Graduated markings to 1000 ml

A lightweight urinal bottle with integrated handle and anatomically shaped opening designed specifically for the female user.

- specifically designed and contoured for women
- opaque bottle, free of sharp edges inside and out
- graduated measurement markings for accurate content measurement
- integrated handle with sturdy grip for easy handling
- durable and easy to clean and disinfect (autoclave at 135°C or sterilise in boiling water or with steam)
- can be used when seated or standing



<b>CODE</b>	BL5421
<b>CAPACITY</b>	1000 ml
<b>DIMENSIONS</b>	155 x 320 x 85 mm
<b>MATERIALS</b>	High Density Polyethylene
<b>WEIGHT</b>	140 g
<b>CONFORMS TO</b>	ARTG 202562
<b>WARRANTY</b>	1 year

DUE TO HYGIENE  
REGULATIONS,  
PRODUCT CANNOT BE  
RETURNED AFTER PURCHASE.