

BetterLiving

Wheeled Walker Push Down

Unique push down
braking mechanism

Height adjustable handles
(810-945 mm)

Lightweight alloy frame with
tube-in-tube construction
for increased strength

Compact fold down design

153 mm (6 inch) castors

Softgrip handles for greater comfort

Foam backrest

Moulded seamless seat with
concealed PVC seat pouch

Removable
basket and bag



Designed for general indoor use, the Push Down Wheeled Walker incorporates a unique brake activation system making it ideal for physically or cognitively impaired users.

The alloy frame makes this wheeled walker especially lightweight, while the tube-in-tube construction significantly increases stability and durability. Six inch castors ensure easy operation and manoeuvrability on interior surfaces and through doorways.

The unique push-down brake mechanism makes the product especially easy to operate for people with physical or cognitive difficulties. Brakes are automatically activated by downward pressure on the handles or by sitting on the seat. Handle height can be adjusted easily to best suit the user, with an extension range of 135 mm to suit users of average weight and height.

The fully moulded PU seat is completely seamless, reducing the potential for skin tear or abrasion, plus foam backrest for added comfort when seated.

BetterLiving

Wheeled Walker Push Down



Cleverly concealed beneath the seat is a PVC pouch, offering additional security for small personal valuables.

With the added convenience of both basket and bag that can be used together (basket fits inside the bag to conceal contents) or individually. Both items are removable, for the security and transport of goods when leaving the wheeled walker unattended.

Easy closing mechanism and compact fold down design for convenient transport or storage.

Colour	Cherry Red
Code	BL1553
Capacity	130 kg
Weight	8 kg
Conforms to	ARTG 176830 ISO 11199-3:2005
Warranty	1 year

HEIGHT	LENGTH	WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR SIZE
945 mm max 810 mm min	610 mm	575 mm	310 mm	360 mm	557 mm	215 mm	153mm 6 inch

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.