Etac Edge, Easy and Smart shower stools

A variety of small and handy shower stools. All stools have soft contoured seats with a non slip pattern and oval drainage holes. These features provide stable and comfortable seating. Choose depending on shower space, function or taste.



Etac Edge shower stool

Fits perfectly in the corner of the shower where the walls provide back stability for the user. Etac Edge can also be used the other way around as a saddle seat which relieves the strain on an aching hip. Also available in a Low version.

Etac Smart shower stool

Just as comfortable whether sitting on the long or short side. Also available with a Swivel pad and in a Low version which can be used as a bath seat. The Low stools are also perfect to use in a care situation.

Etac Easy shower stool

Easy to use, easy to move. A round seat is always in the right position. Also available with a swivel pad.



Features



Easy to fit

The push button lock function makes it easy to both fit and remove the legs.



Easy and safe The oval shape prevents the telescopic legs from turning when adjusting the height.



Stability

To ensure stability, the legs are angled outwards. This means that the higher the seat, the larger the footprint.



Ferrules with a firm grip

The ferrules squeeze away the water and adapt perfectly to the floor.

Accessories



Seat pads Edge: 81801132

 \checkmark

Smart: 80209502 Easy: 80209501 The non-slip seat pads provide an instant warm feeling and are water resistant.



Etac Easy and Smart with Swivel pad enables easy turning without any friction. Removable for cleaning. Material: Polyethylene. Cleaning at max. 60°C (140°F).

	Etac shower stools	Grey		rey, with wivel pad	Weight Edge: 2.6 kg (5.7 lbs)	
<u> </u>	Edge shower stool	818010	81801040		Edge Low: 2.4 kg (4.4 lbs) Easy: 2.8 kg (6.2 lbs)	
∕⊻	Easy shower stool	819010	8	1901050	Smart: 3.2 kg (7.1 lbs)	<
\checkmark	Smart shower stool	819510	81951030 8 1	1951050	Smart Low: 3 kg (6.6 lbs)	
	Dimensions Edge Total width: 52 cm (20½") Seat width: base max 45 cm Seat height Edge: 42–57 cm (Seat height Edge Low: 34.5–42 cm (13½–16½")	. ,	Dimensions Smart ootprint: 45.5–46,5 x 41–4. 18–18¼ x 16¼–17") eat width: 43–38 cm (17–1 eat height: 42–57 cm (16½ eat height Smart Low: :4.5–42 cm (13½-16½")	5")	Colours Grey Max. user weight Smart and Easy: 150 kg (330 lbs) Edge: 130 kg (287 lbs)	
	Dimensions Easy Footprint: 49,5–53,5 cm (19) Seat width: Ø40 cm (Ø15¾") Seat height: 42–57 cm (16½)	Aaterials eat: Polypropylene egs: Aluminium, polyamide errules: TPE wivel pad: Polyethylene	or with a	product with a solvent-free cleaning agent (pH level 5–9), 70% disinfectant solution. Can be decontaminated at max. 5°F if so required.	