

# Bassett Lumbar Support

Sturdy steel frame

More positioning clips  
on Heavy Duty model



Contoured, open woven mesh  
for airflow and support

Velcro strap for secure  
seat attachment



A lightweight back rest designed for use in cars, chairs and wheelchairs. It features a curved metal frame with a contoured open mesh to support the back.

Specifically recommended for people with arthritic conditions, providing excellent relief for back pain, neck pain, sciatica and posture problems.

The 'D' shaped metal frame rests firmly and securely against the chair, and supports correct lumbar posture. The mesh promotes good air circulation to prevent undue perspiration, while providing maximum support to the lower back. Extremely light and portable for easy transport and storage. Includes a strap for secure seat attachment.

The Heavy Duty model includes a more pronounced lumbar support with more closely woven cord and a greater number of plastic clips.

An optional cover made from polar fleece is also available, helpful to stop back pocket buttons from catching in the mesh.

ITEM	Standard	Heavy Duty	Cover
CODE	Standard-Support	Heavy-Duty-Support	Cover-for-Bassett
HEIGHT	440 mm	440 mm	
WIDTH	390 mm	390 mm	
DEPTH	110 mm	135 mm	
MATERIALS	Powder coated mild steel frame Polyester knotless net mesh and webbing, velcro and nylon buckle		
WEIGHT	850 g	950 g	
WARRANTY	1 year		

