

Etac Relax shower seat

Scandinavian design for the modern bathroom.

Attention to detail

Etac Relax has been designed to ensure a secure and enjoyable shower experience. Design, function and safety are extremely important to us, and this is presented here in the form of a discreet, wallmounted shower seat. The shape of the seat offers the option of varied sitting positions, resulting in improved access. It is always close at hand, but as it can be folded upwards it never takes up any unnecessary space.

Adapted to your needs

Etac Relax comes in two colours and six different models, which have all been designed to create options specifically for you.





	Etac Relax shower seat	
	- White	81703010
	- White, with arm supports	81703020
$\langle \nabla \rangle$	- White, with supporting legs	81703030
	- White, with arm supports and back support	81703060
$\overline{m{ee}}$	- White, with arm supports and supporting legs	81703040
	- White, with arm supports, back support and	81703050
	supporting legs	
	- Volcano grey	81708000
	- Volcano grey with arm supports	81708010
	- Volcano grey with supporting legs	81708020
	- Volcano grey with arm supports and	81708050
	back support	
	- Volcano grey with arm supports and	81708030
	supporting legs	
	- Volcano grey with arm supports,	81708040
	back support and supporting legs	

Colour

White, Volcano grey

Maximum user weight

With/without supporting legs: 150/125 kg

Material

Seat; polypropylene Arm supports; polypropylene and TPE Legs; aluminium and stainless steel, Ferrules; TPE



 ϵ

Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max $85^{\circ}C$ / $185^{\circ}F$ if so required.

Features



Large, stable seat Large and stable seat when folded down and extremely discreet when folded up.




Easy-grip arm support The shape and material of the folding arm supports provide a safe and secure grip.



Functional design The shape of the seat offers the option of varied sitting positions, resulting in improved access.



Easy to fit As there are only three fixing points, installation is both quick and easy.



Screws/fastener are not included. May be mounted on walls with sufficient load-bearing capacity only. Please ensure that suitable fixtures are used. Mounting of the seat must be carried out by a qualified person. More information in Etac Relax user manual at www.etac.com.

Accessories



Back support, soft

Soft heat reflective back support. Attached to the wall with velcro. 81704060 White



Supporting legs

For extra stability and user weight up to 150 kg (330 lbs). Height adjustable 45-55 cm (17¾"-21½").

White Volcano grey





For mounting on wood/plywood or concrete walls. The kit includes 3 screws and plugs

84006130





81704050

81708070