

This convenient, portable blood pressure monitor wraps around the wrist and uses the oscillometric principle to simply and accurately measure systolic pressure, diastolic blood pressure and pulse rate.

With an easy fitting wrist cuff and one touch operation, the unit automatically inflates to the correct pressure, using oscillometric technology to simultaneously measure the arterial blood pressure and pulse rate.

A 90 memory storage capacity allows the review of the last 90 readings.

Also features low battery indicator and automatic power down after 3 minutes of no operation. Requires 2 x AAA batteries (included) and suits wrist circumference of 125 to 205 mm.

The compact unit is lightweight and easy to carry when travelling or for daily use.

CODE	BSSDS733	
CUFF	125 > 205 mm circumferenc	e
HEIGHT	25 mm	
WIDTH	61 mm	
DEPTH	62 mm	
WEIGHT	190 g with batteries	
POWER	2 x AAA 1.5V batteries	\ /
ACCURACY	Blood Pressure ±3 mmHg	
	Pulse Rate ±5%	_
CONFORMS TO	ARTG165748	\
WARRANTY	1 year	

FOR PEOPLE WITH HIGH BLOOD PRESSURE. CONSULT YOUR MEDICAL PROFESSIONAL TO EVALUATE THE READINGS, CHECK YOUR DEVICE PERIODICALLY FOR ACCURACY.