

REGAL S344A





FEATURES & BENEFITS

- Loop handlebars offer comfortable finger or thumb throttle operation
- 4 wheel adjustable suspension system allowing for more rider control and comfort
- Large 14" wheels which provide greater ground clearance
- Larger 4 pole motor and larger controller delivering much more power for climbing steep driveways or hills
- Built-in half speed safety switch which automatically cuts in when turning to provide extra safety
- New improved, fully adjustable leather-feel seat, allowing maximum rider comfort
- Separate emergency braking system for added security

Speed	10 -15 kph (programmable)
Range	30-50 km
Motor	650W / 2200W, 4 Pole
Controller	Rhino 160 amp
Battery	12V / 55-75 a/h x2 pieces
Charger	8 amp off-board
Gradient	10-12°
Front Wheels	14" pneumatic
Rear Wheels	14" pneumatic
Weight capacity	227 kg
Brake	Intelligent, regenerative, electromagnetic brakes
Length	152 cm
Width	71 cm