

Shoprider

Rocky 6

ASK FOR THE BRAND SHOPRIDER



ROCKY 6
Fat Tyre

Available Colours



SHOPRIDER®
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

Rocky 6

SHOPRIDER
www.shoprider.com.au
mobility

The SHOPRIDER™ **Rocky 6 Fat Tyre** is the flagship of the SHOPRIDER range, it redefines the concept of the Luxury in mobility scooters. The **Rocky 6 Fat Tyre** comes standard with 330mm wheels, Double Pillow Comfort Seat and all the features you expect from the SHOPRIDER range, fully rotating lockable seat, comfortable suspension (adjustable on the rear), gas adjustable steering tiller and multi position fold down arm rests. The reliability of this impressive scooter by SHOPRIDER is further enhanced by the strong 4 pole motor providing rear wheel and electromagnetic brake.

SPECIFICATIONS	ROCKY 6 Fat Tyre
Length	1430mm
Width	680mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size - 4 Pole	2.0Hp
Maximum Speed	10Kph
Battery capacity min	2 x 60 Amp
Maximum range#	up to 50km
Wheel Size	330mm
Turning Radius	1800 mm
Weight (with battery)	146kg

FEATURES:


- Double Pillow Comfort Seat
- Delta Steering Tiller
- Gas Adjustable Steering Tiller
- Adjustable width arm rests
- Adjustable rear suspension
- Sliding Rotating Seat
- Variable Seat Height
- Battery Power indicator
- Double Head Lights
- Rear Lights & Indicators
- Hazard Lights
- Front & Rear Bump Pads
- High & Low Speed
- Shopping Basket
- Comprehensive Warranty**

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

SHOPRIDER
www.shoprider.com.au
mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au