

ALUMINIUM WHEELCHAIR RANGE





Travelite Transport Chair



Aluminium Travel Chair



Aluminium Travel Chair Plus



Travel Chair in a Bag



Lightweight Aluminium



Ultra Lightweight Aluminium



XS



XS₂



SD₂





TRAVELITE TRANSPORT CHAIR

Sturdy yet lightweight, this occasional use chair is ideal for many outings from shopping trips to holidays.



- Folds to approx 1/3 the size of a typical folded wheelchair.
- · Ideal for holidays, days out and shopping trips.
- · Flip up padded desk arms.
- Multiple steel cross braces for added strength.
- Aluminium frame in gun metal grey.
- 8" puncture proof front and rear tyres.
- · Lightweight travel bag with wheels.
- 115kg (18st) weight capacity





ALUMINIUM TRAVEL CHAIR

The budget travel chair is lightweight and folds to a compact size making it ideal for days out etc.



"A generous 19" seat width."

- · Lightweight aluminium frame.
- A half-folding back for easy transport and storage.
- · Comfortable padded armrests.
- A 19" width seat canvas in nylon with lift up straps and lap belt.
- · Parking brakes.
- · 8" puncture proof tyres.
- Quick release swing out footrests.
- · Stepper tube.
- · Weighs 19lb (8.6kg) including all components.
- Cranked push handles to reduce storage space.
- 115kg (18st) weight capacity

ALUMINIUM TRAVEL CHAIR PLUS

This upmarket travel chair is ideal for the user who requires additional features.



"Easier to control thanks to cable brakes."

- · Padded back canvas with phone and map pockets.
- · A half-folding back for easy transport and storage.
- · Comfortable padded armrests.
- A 19" width seat canvas in nylon with lift up straps and lap belt.
- · Parking brakes.
- 8" front and 12" rear puncture proof tyres.
- · Quick release swing out footrests.
- Flip up stepper tube.
- Weighs 25lb (11.5kg) including all components.
- Compact design with cranked push handles to reduce storage space.
- 115kg (18st) weight capacity

TRAVEL CHAIR IN A BAG

The Enigma Travel Chair is lightweight at only 19lb making it ideal for transportation.



"Comes complete with bag as standard."

- · Quick release / swing out height adjustable footrests.
- · Half folding back mechanism.
- 19" width seat canvas in nylon upholstery with lift-up straps and lap belt.
- · Parking brakes.
- 6" front and 12" rear puncture proof tyres.
- · Flip up stepper tube.
- · Desk friendly compact armrests.
- · Padded back canvas with phone and map pockets.
- · Comes complete with a storage bag and is ideal for transportation.
- · Metallic silver finish.
- · Tool kit included.
- · Total weight only 8.6kg (19lb).
- Weight without footrests only 7.4kg (16.3lb).
- 115kg (18st) weight capacity

LIGHTWEIGHT ALUMINIUM

The Lightweight Aluminium Wheelchair offers a value for money introduction to the Enigma Aluminium Wheelchair range without compromising on quality.



"Our most popular aluminium wheelchair."

- 18" seat width with a nylon seat & back canvas.
- 8" Puncture proof front tyres.
- · Half folding back mechanism.
- Full length armrests with nylon upholstered arm pads.
- · Quick release, swing away, height adjustable footrests.
- · Folds up easily for storage and transportation.
- · Weighs just 8.5kg (18.7lb) in carry format.
- 115kg (18st) weight capacity.

SELF PROPEL MODEL

- · Quick release aluminium 24" wheels with hand rim and grey pneumatic tyres.
- · Weighs only 14.5kg (32lb) in total.

TRANSIT MODEL

- Weighs only 12.4kg (27.3lb) in total.
- Includes cable brakes, which an attendant may use to reduce speed.
- 12.5" solid puncture proof rear tyres.

ULTRA LIGHTWEIGHT ALUMINIUM

The Ultra Lightweight Aluminium chair was developed to include all the features of the LAW range with an additional focus on design.



"This ultra light chair features a stylish new design."

- Silver aluminium frame with single cross brace.
- · Grey two-tone upholstery with rear pocket.
- · Swing out quick release footrests with height adjustment.
- · Half folding back mechanism.
- · Desk style armrests with fabric side panels and pocket.
- Two stepper tubes.
- Plastic front forks with 20cm (8") x 3.2cm (1") grey PU tyres.
- 115kg (18st) weight capacity.

SELF PROPEL MODEL

- · Long reach brakes.
- · Spoked rear wheels with 24" x 13/8" grey solid puncture proof tyres.

TRANSIT MODEL

- · Cranked push handles for compact folding.
- · Long reach brakes and rear cable brakes.

XS ALUMINIUM

The XS Aluminium was developed to provide more features as standard with added comfort.



"Due to popular demand now available in a 20" seat width and in Spirit Red."

- Choice of 24" quick release wheels (self-propel) or 12.5" fixed wheel (transit).
- Padded tension adjustable back and seat canvasses for comfort.
- · Comes with deluxe seat and back cushions.
- · Quick release swing in / out footrest with full height and footplate angle adjustment.
- · Flip back and removable armrests.
- · Half folding back mechanism.
- 18" or 20" seat width x 17" depth seat canvas.
- · Height adjustable anti-tip wheels.
- · Adjustable rear wheel position.
- · Weighs under 9kg in carry format.
- · Crash tested.
- · Solid puncture proof tyres.
- · Available in blue and sprit red.
- 115kg (18st) weight capacity.

Colour Options





XS₂

Our new and improved XS Wheelchair uses a similar frame to the XSAWC Wheelchair but adds many new features. The chair is available in a vivid Silver in 16", 18" and 20" seat widths.



"We took one of our most popular chairs and simply made it better"

- One touch lap belt as standard on both self propel and transit versions.
- · Solid puncture proof front and rear tyres.
- · Length adjustable armrests.
- · Long reach brakes on self propel.
- · Improved footrest release catch.
- Choice of 24" quick release wheels (self-propel) or 12.5" fixed wheel (transit).
- · Padded tension adjustable back canvas for comfort.
- · Comes with seat and back cushions.
- Quick release swing in / out footrest with full height and footplate angle adjustment.
- · Half folding back mechanism.
- 16", 18" or 20" seat width x 17" depth seat canvas.
- · Height adjustable anti-tip wheels.
- · Adjustable rear wheel position.
- · Crash Tested.
- · Weighs under 9kg in carry format.
- 135kg (21st) weight capacity.

SD₂

Introducing the NEW SD2 Wheelchair, which is an upgrade in design and features and will replace our old model, the Super Deluxe chair.



"For those who require maximum features."

- · Elegant metallic black aluminium frame.
- Long reach brakes.
- Flip up and removable armrests.
- · Push button arm pad adjustment height and length.
- · Improved swing in/out footrests.
- Half fold back mechanism.
- · Tension adjustable back with cushion.
- · Stepper tube with integrated anti tip wheels.
- Attendant cable brakes on Transit models.
- · Reflectors included.
- · Quick release rear wheels on Self Propel models.
- · Black pneumatic tyres.
- · One touch lap belt.
- Available in 16", 18", 20" and 22" seat widths.
- Storage pocket on back of the chair.
- · Comes with seat and back cushions.
- 135kg (21st) weight capacity.

SPIRIT WHEELCHAIR

The Spirit wheelchair uses the latest in wheelchair design to produce an elegant yet highly practical wheelchair. The sculptured frame is constructed from lightweight aluminium for durability and responsiveness.



The ultimate in style and design'

- · One-touch adjustable armrests for easy transfer and comfort.
- · Folding frame with removable footrests and half-folding back mechanism.
- · Padded seat and back cushion.
- · A tension adjustable back canvas.
- · Puncture-proof solid front and rear tyres.
- · Long reach brake handles.
- · Height adjustable anti-tip wheels.
- · Quick release rear wheels.
- · Available in electric blue, spirit red and silver.
- · Available with MAG wheels (red only).
- 115kg (18st) weight capacity.

Colour Options







K - CHAIR

The Enigma K Chair has been developed to provide a new level of comfort for wheelchair users. It features adjustable suspension that can be tailored for the user and the environment thus reducing the effects of any uneven floor surface.



"Adjustable suspension for optimum comfort."

- · Half folding back for easy storage.
- Padded seat and back cushion.
- · Tension adjustable back canvas.
- · Rear shock absorbers, adjustable for ride and comfort.
- Shock absorbing front castors with 7" puncture proof tyres.
- 24" rear wheels (quick release) with hand rims, black puncture proof solid tyres and 6-position axle.
- One touch height adjustable armrests.
- · Height adjustable anti-tip wheels.
- · Swing in / swing out footrests with new easy to use pull up handle.
- Stylish ergonomic hand grips.
- · Long reach brakes.
- · Available in midnight blue or spirit red.
- 115kg (18st) weight capacity.

Colour Options





SPECIFICATION

CODE	DESCRIPTION	WIDTH (OPEN)	WIDTH (CLOSED)	LENGTH	НЕІСНТ	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
TC005	Travelite Transport Chair 18"	57cm (22")	29cm (12")	84cm (33")	95cm (37")	45cm (18")	45cm (18")	48cm (19")	10kg (21lb)	115kg (18st)
ATC19-BK	Aluminium Travel Chair 19"	56cm (22")	23cm (9")	84cm (33")	99cm (39")	48cm (19")	41cm (16")	48cm (19")	8.6kg (19lb)	115kg (18st)
TC002SIL	Lightweight Travel Chair In a Bag 19"	56cm (22")	23cm (9")	87cm (34")	90cm (35")	48cm (19")	39cm (15")	46cm (18")	8.6kg (19lb)	115kg (18st)
EXP19BL	Aluminium Travel Chair Plus 19"	56cm (22")	31cm (12")	89cm (35")	99cm (39")	48cm (19")	40cm (16")	48cm (19")	11.5kg (25lb)	135kg (21st)
LAWCoo1	Lightweight Aluminium Wheelchair 18" Self Propel	61cm (24")	24cm (9")	106cm (42")	90cm (35")	46cm (18")	40cm (16")	48cm (19")	15kg (32lb)	115kg (18st)
LAWC002	Lightweight Aluminium Wheelchair 18" Transit	55cm (22")	29cm (11")	96cm (38")	91cm (36")	46cm (18")	40cm (16")	48cm (19")	12kg (27lb)	115kg (18st)
LAWC007A	Ultra Lightweight Aluminium 18" Self Propel	62cm (24")	32cm (13")	107cm (42")	93cm (37")	42cm (17")	39cm (16")	52cm (20")	13kg (29lb)	115kg (18st)
LAWC008A	Ultra Lightweight Aluminium 18" Transit	61cm (24")	28cm (11")	95cm (38")	97cm (36")	42cm (17")	39cm (16")	52cm (20")	12kg (26lb)	115kg (18st)
LAWConA	Ultra Lightweight Aluminium 20" Self-Propel	71cm (28")	32cm (12.6")	109cm (43")	93cm (37")	51cm (20")	39cm (16")	52cm (20.5")	13kg (29lb)	125kg (20st)
LAWC012A	Ultra Lightweight Aluminium 20" Transit	68cm (27")	28cm (11")	103cm (41")	97cm (38")	51cm (20")	39cm (16")	52cm (20.5")	n.6kg(25.6lb)	125kg (20st)
XSAWCTS18ST	XS Aluminium 18" Transit	63cm (25")	26cm (10")	92cm (36")	93cm (36·5")	45cm (18")	43cm (17")	52cm (20·5")	15kg (32lb)	115kg (18st)
XSAWCSP18ST	XS Aluminium 18" Self Propel	63cm (25")	26cm (10")	103cm (40·5")	93cm (36·5")	45cm (18")	43cm (17")	52cm (20·5")	17kg (37lb)	115kg (18st)
XSAWCTS20ST	XS Aluminium 20" Transit	68cm (27")	26cm (10")	92cm (36")	93cm (36·5")	50cm (20")	43cm (17")	52cm (20·5")	15kg (32lb)	115kg (18st)
XSAWCSP2oST	XS Aluminium 20" Self Propel	68cm (27")	26cm (10")	103cm (40·5")	93cm (36·5")	50cm (20")	43cm (17")	52cm (20·5")	17kg (37lb)	115kg (18st)

XS2SP16SIL	XS2 Aluminium 16" Self-Propel	59cm (23")	26cm (10")	103cm (40.5")	93cm (37")	40cm (16")	43cm (17")	52cm (20.5")	17kg (37lb)	135kg (21st)
XS2TS16SIL	XS2 Aluminium 16" Transit	59cm (23")	26cm (10")	92cm (36")	93cm (37")	40cm (16")	43cm (17")	52cm (20.5")	17kg (37lb)	135kg (21st)
XS2SP18SIL	XS2 Aluminium 18" Self-Propel	63cm (25")	26cm (10")	103cm (40.5")	93cm (37")	45cm (18")	43cm (17")	52cm (20.5")	17kg (37lb)	135kg (21st)
XS2TS18SIL	XS2 Aluminium 18" Transit	63cm (25")	26cm (10")	92cm (36")	93cm (37")	45cm (18")	43cm (17")	52cm (20.5")	15kg (32lb)	135kg (21st)
XS2SP2oSIL	XS2 Aluminium 20" Self-Propel	68cm (27")	26cm (10")	103cm (40.5")	93cm (37")	50cm (20")	43cm (17")	52cm (20.5")	17kg (37lb)	135kg (21st)
XS2TS2oSIL	XS2 Aluminium 20" Transit	68cm (27")	26cm (10")	92cm (36")	93cm (37")	50cm (20")	43cm (17")	52cm (20.5")	15kg (32lb)	135kg (21st)
SD2SP/TS16BLK	Super Deluxe 216" Self Propel/Transit Black	68cm (27")	28cm (11")	103cm (40.5")	93cm (37")	40cm (16")	41cm (16")	57cm (22")	18kg (39lb)	135kg (21st)
SD ₂ SP/TS ₁₈ BLK	Super Deluxe 218" Self Propel/Transit Black	68cm (27")	28cm (11")	103cm (40.5")	93cm (37")	45cm (18")	41cm (16")	57cm (22")	18kg (39lb)	135kg (21st)
SD2SP/TS20BLK	Super Deluxe 2 20" Self Propel/Transit Black	68cm (27")	28cm (11")	103cm (40.5")	93cm (37")	52cm (20")	41cm (16")	57cm (22")	19kg (41lb)	135kg (21st)
SD2SP/TS22BLK	Super Deluxe 2 22" Self Propel/Transit Black	68cm (27")	28cm (11")	103cm (40.5")	93cm (37")	56cm (22")	41cm (16")	57cm (22")	19kg (41lb)	135kg (21st)
XSES18RD	Spirit Wheelchair 18" - Spirit Red	71cm (28")	32cm (12·5")	97cm (38")	91cm (36")	45cm (18")	40cm(15·5")	48cm (19")	15kg (32lb)	115kg (18st)
XSES18BL	Spirit Wheelchair 18" - Electric Blue	71cm (28")	32cm (12·5")	97cm (38")	91cm (36")	45cm (18")	40cm(15·5")	48cm (19")	15kg (32lb)	115kg (18st)
XSES18SI	Spirit Wheelchair 18" - Silver	71cm (28")	32cm (12·5")	97cm (38")	91cm (36")	45cm (18")	40cm(15·5")	48cm (19")	15kg (32lb)	115kg (18st)
XSESMAG18RD	Spirit Wheelchair 18" with MAG Wheels - Red	71cm (28")	32cm (12·5")	97cm (38")	91cm (36")	45cm (18")	40cm(15·5")	48cm (19")	15kg (32lb)	115kg (18st)
XSSW18RD	Enigma K-Chair 18" - Spirit Red	64cm (25")	28cm (11")	108cm (42·5")	97cm (38")	45cm (18")	43cm (17")	49cm (19")	20kg (44lb)	115kg (18st)
XSSW18BL	Enigma K-Chair 18" - Midnight Blue	64cm (25")	28cm (11")	108cm (42·5")	97cm (38")	45cm (18")	43cm (17")	49cm (19")	20kg (44lb)	115kg (18st)

Comfort, Style, Performance & Reliability.

We understand the true value of innovation within this market place; and know that our Wheelchair development will have a direct influence on improving an individual's life.



All sizes and weights stated are nominal. All errors and omissions excepted.

Due to the limitations of photography and the printing process colours shown may not be 100% accurate.