

NEW



Yes, you can.

## Invacare® **Comet™ Heavy Duty**



### Well-proportioned, Comfortable, Supporting: The *Invacare Comet HD*

The Invacare **Comet HD** is part of the Invacare Comet family, renowned for its safe, fast and enjoyable drive. This rugged outdoor vehicle is ideal for users weighing up to 220 kg. Packed with safety features, this powerful scooter will take you to your destination in the most direct and secure way. The large contoured extra padded seat together with four individual adjustable suspension units ensures an exceptional driving experience. The **Comet HD** has been designed to allow for additional weight capacity and has passed thorough tests in the German independent TÜV test laboratory.



## ➤ Durable Design and Power in One


Even for users with a weight of up to 220 kg, safety is never a concern with the **Comet HD**. All necessary precautions for guaranteeing the stability and safety of the scooter have been taken into consideration during the design process.

## ➤ That Comfortable Feeling

The XL contoured seat has been developed in conjunction with therapists to ensure optimum seating comfort. The seat base and height- and width-adjustable armrests have been reinforced to allow for additional stresses and can easily be used as supports when getting on and off the scooter. Seating comfort can be increased by using the pneumatically adjustable lumbar support.



## Features



**Reinforced seat column**  
It ensures stability and safety while seated.



**Reinforced armrest**  
It allows for the additional stress of getting on and off the **Comet HD**.



**Unique suspension system**  
It has proven itself in the professional mountain biking industry and allows for simple independent user adjustment.

## Technical Data



66cm



53cm



44-51.5cm



61cm



66cm



147.5cm




148kg\*



220kg



(2) 12V 75Amp



326cm



10cm



13.5cm



10°



Cont: 650W  
Peak: 1800W



50km



10 km/h

Model:	4 Wheel
Colours:	Blue, Silver & Red
Wheels:	Front & Rear 31 cm/13" Pneumatic
Suspension:	Front & Rear Spring & Adjustable (Heavy Duty)
Seat Type:	XL contoured seat, Reinforced seat column & armrest, Extra padding, ergonomic, recliner lever, & headrest
Std. Batt. Charger:	24V / 8 Amp (off-board)
Electronics:	160 Amp Rhino Controller

Compliance: TGA-registered (ARTG# I 19338) CE Marked TUV certified – tested to European Standard (DIN EN 12184:1999)

\*Without batteries

## Your Local Dealer:



Yes, you can.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Invacare Australia Pty Ltd ABN 45 074 676 378.  
PO Box 5002, 1 Lenton Place, North Rocks, NSW 2151, Australia.  
Freephone: 1 800 460 460, Fax: 02 8839 5353,  
Email: sales@invacare.com.au, www.invacare.com.au

Invacare New Zealand,  
PO Box 62-124, 4 Westfield Place, Mt Wellington, Auckland, New Zealand.  
Freephone: 0508 468 222, Freefax: 0508 807 788  
Email: sales@invacare.co.nz, www.invacare.co.nz

Invacare Asia Ltd,  
PO Box 5002, 1 Lenton Place, North Rocks, NSW 2151, Australia.  
Phone: +61 2 8839 5330, Fax +61 2 8839 5331,  
Email: asiasales@invacare.com, www.invacare.com