

## Folding Walking Frames



### Folding Walking Frames

Model RG4070WB -	32"- 38"	Height Adjustable 5" Wheel with Rubber Grip
Model RG4070MD -	32"- 38"	Height Adjustable, Rubber Grip
Model RG4070MA -	29.5" – 33.5"	Height Adjustable, Foam Grip

### Dimensions

- Max user capacity 150kg
- Walker net weight 3.0kg – 3.2kg

### Features

- Aluminium anodized frame
- Multi height adjustment on each model
- Soft foam or rubber handgrips

#### Optional Extra

- Swivel Wheels
- Wheels & Glides

---

Phone: 08 9248 4180  
Fax: 08 9248 9452  
E-mail: [orders@redgumbrand.com.au](mailto:orders@redgumbrand.com.au)