

Montreal Wall Hugger Lift Chair

The Montreal Wall Hugger Lift Chair has been designed to offer superior seating comfort as it is fitted with a pocket spring seat base plus the added functionality of the space saving wall hugger operation. The wall hugger operation has a legs up only action and also allows the chair to be positioned 100mm from the wall with full operation available on the chair which is ideal for smaller rooms. Each chair is fitted with quality actuators and are available in either Macro Suede or Faux Leather covering.



Specifications

Maximum User Weight	150kg
Overall Width	845mm
Overall Depth Upright	820mm
Overall Height	1145mm
Seat Width	490mm
Seat Depth	520mm
Seat Height	500mm
Back Rest Height	695mm
Arm Rest Height	160mm

Colour Availability

LC0501	Chocolate Brown	Macro Suede
LC0502	Beige	Macro Suede
LC0503	Dark Blue	Macro Suede
LC0504	Brown	Faux Leather
LC0505	Taupe	Faux Leather
LC0506	Black	Faux Leather

Warranty Information

Steel Base	3 Years
All other components	1 Year

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au