

REDGUM BRAND
Daily Living and
Lifestyle Products

European Designed Deluxe Domestic Care Electric Bed

RGBDK5 LUNA



Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au

REDGUM BRAND
Daily Living and
Lifestyle Products

*REDGUM Brand is associated with Shoprider Mobility
www.redgumbrand.com.au or www.shoprider.com.au Ask for the Brand SHOPRIDER*

European Designed Deluxe Domestic Care Electric Bed

RGBDK5 LUNA



Features

- European Quality Design
- 4 Way Sectional Laying Surface Adjustment – Height, Back, Foot and Thigh
- Anti Trendelenburg
- Strong Metal Slatted Base
- Ergonomically Designed Hand Controller with Large Control Buttons.
- Lock Out feature on each function of the Hand Controller.
- Integrated Wooden Twist Proof Side Rails.
- Up-Lifting Pole Fitted as standard
- Locking Castor Wheels.
- Delivered with Transport Kit for Easy Movement.
- Quality Linak Actuators

Dimensions

Mattress Surface	90 x 200cm
Overall Size	103 x 215cm
Max User Weight	140kg
Net Weight of Bed	80kg
Safe Working load	175kg
Castor Wheel Size	100mm
Height Adjustment	40 x 80cm