

Foldable Shower Chair

RG9400



The REDGUM Brand Foldable Shower Chair has been designed to easily fold for storage when not required. This foldable shower chair has tough aluminium painted frame with stainless steel hinge pins with a polyurethane self draining seat. The legs are easily adjusted to suit

multiple seat height requirements and can be adjusted to accommodate uneven floors. Each chair leg is fitted with durable rubber stoppers to minimise any movement. The foldable shower chair is finished a pearlescent white finish.

Dimensions

- Seat width 470mm
- Leg height adjustment range 460mm to 530mm
- Max user capacity 100kg

Features

- Easy fold design
- Self draining polyurethane seat
- Robust aluminium construction with stainless steel fitting

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au