

Aluminium Shower Stool

with Optional Backrest
Australian Standards Approved



Features

- Australian Standards Approved
- Lightweight Aluminium Construction
- Anodised Frame
- Polyurethane Seat with Drain Points
- Non Slip Seat Finish
- Multi Slot Legs for Uneven Floors
- Soft Padded Arm Rests
- High Grade Rubber Feet
- Rust Free Fittings

Dimensions

- Seat Width -455mm
- Seat Depth -330mm
- Seat Height - 470mm -595mm
- Max User Capacity - 105kg
- Stool Weight - 2.8kg



Options Available

RG6002

- Removable Backrest
- Multi Position
- Easy Clean Polyurethane Material
- Lifting Handle in Backrest

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au