

Heavy Duty Shower Stool

RG9012 / RG9012R



The REDGUM Brand Round Bath Stool is available in two models, one has a unique rotating padded seat and the other has a fixed padded seat. These stools have been designed to withstand the rigours of the bathroom environment. The stool is extremely stable with strong aluminium legs that can be adjusted to suit multiple heights, each leg also has large rubber stoppers for added stability

Dimensions

- Seat width 325mm
- Leg height adjustment range 480mm to 670mm
- Max user capacity 150kg

Features

- Strong rust free construction
- Removable legs for easy storage
- Rotating seat on RG9012R

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au