

REDGUM BRAND
Daily Living and
Lifestyle Products

European Designed Standard Domestic Care Electric Bed

RGBDL5



Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au

REDGUM BRAND
Daily Living and
Lifestyle Products

*REDGUM Brand is associated with Shoprider Mobility
www.redgumbrand.com.au or www.shoprider.com.au Ask for the Brand SHOPRIDER*

European Designed Standard Domestic Care Electric Bed

RGBDL5



Stylish Wooden & Metal
Highlighted Finish



Ratchet Height Adjustment
to Lower leg



Uplifting
Pole fitted
as standard

Features

- European Quality Design & Manufacture
- 4 Way Sectional Laying Surface Adjustment – Height, Back, Foot and Thigh
- Strong Metal Mesh Base
- Ergonomically Designed Hand Controller with Large Control Buttons.
- Safety Lock Feature on Hand Controller .
- Energy Saving Mode allows the Power to the Bed to go to Standby When Not in Use.
- Integrated Wooden Twist Proof Side Rails.
- Up-Lifting Pole Fitted as standard
- Locking Castor Wheels

Dimensions

Mattress Surface	90 x 200cm
Overall Size	105 x 219cm
Max User Weight	150kg
Net Weight of Bed	90kg
Safe Working load	185kg
Castor Wheel Size	100mm
Height Adjustment	40 x 80cm

Options

- Foam Mattress with Cotton Cover
- Ribbed Comfort Mattress
- Moisture Resistant Mattress Cover
- Bed Lengthening Kit by 200mm