

# Afiscooter C4

The outstanding design of the Afiscooter-C is highlighted by these user comfort features:

fully adjustable swiveling orthopedic seat that gives easy access and an enjoyable ride, as well as the DELTA tiller that provides easy steering and is fully adjustable to accommodate any driver.

A comfortable drive is guaranteed with the advanced full suspension system.

In addition, the Afiscooter-C features a large rear basket for convenient storage.

The four wheel model allows for a smooth and comfortable ride along with improved stability, indoors and outdoors.

Available Colours



## **Afikim Electric Vehicles**

has been developing and manufacturing advanced, electrically-powered mobility scooters since 1978!

Afiscooters' mobility scooters are designed with reliability and driver safety in mind first and foremost alongside beautiful, and thoughtful design.

Afikim Electric Vehicles' worldwide distribution network offers year-round local technical support and meets European ISO 9001-2000, EN 12184 certification, as well as US safety standards with FDA approval.

**RIVERINA**  
**SCOOTERS**  
& Mobility



8 Blake Street, Wagga Wagga NSW 2650

Phone 02 6921 4444

[www.riverinascooters.com.au](http://www.riverinascooters.com.au)

## Features

- ▶ A fully adjustable orthopedic seat provides comfort while driving and swivels for ease of boarding
- ▶ Fits standard car lifts for the traveling user
- ▶ An optional canopy is available for protection from the sun
- ▶ Easy tiller adjustment for any driver



All operational controls including lights, horn and battery display are conveniently located in the center of the one-hand operation Delta handlebar



Plenty of storage space in the standard: front basket and optional large rear basket



## Specifications

|                                    |                   | Afiscooter-C4                       |         |
|------------------------------------|-------------------|-------------------------------------|---------|
| <b>Overall length</b>              |                   | 52"                                 | 132 cm  |
| <b>Overall width</b>               |                   | 26"                                 | 67 cm   |
| <b>Overall height</b>              | Without canopy    | 50"                                 | 126 cm  |
| <b>Ground clearance</b>            |                   | 4"                                  | 10 cm   |
| <b>Weight</b>                      | With batteries    | 269 lb                              | 122 kg  |
|                                    | Without batteries | 205 lb                              | 93 kg   |
| <b>Travel range</b>                |                   | 25 m                                | 40 km   |
| <b>Weight capacity Inc. driver</b> |                   | 330 lb*                             | 150 kg  |
| <b>Maximum speed</b>               |                   | 9.3 mph                             | 15 Km/h |
| <b>Maximum climbing angle</b>      |                   | 16 %                                |         |
| <b>Turning radius</b>              |                   | 67"                                 | 170 cm  |
| <b>Seating</b>                     |                   | 18"                                 | 46 cm   |
| <b>Tires type</b>                  |                   | Pneumatic                           |         |
| <b>Tires size</b>                  | Front             | 4.1/3.5X6"                          |         |
|                                    | Rear              | 4.1/3.5X6"                          |         |
| <b>Battery type</b>                |                   | Two 12V 45-50A/h                    |         |
| <b>Power unit</b>                  |                   | Rear wheel drive, 24V DC motor      |         |
| <b>Color</b>                       |                   | Metallic Silver / Blue / Cherry Red |         |

Tested by manufacturer, based on actual parameters.

\* US Model


**AFIKIM** | 1-800-809-3010 (US only)  
 ELECTRIC VEHICLES | www.afiscooters.com

