

## Trakka Wheelchair Wide Wheels

RG18TW & RG20TW



## Features

- 45mm Wide tyres / Solid Filled
- 23" Aluminum wheels & Push rims
- Strong steel construction
- Powder coated folding frame
- Padded arm rests
- Lift up armrests
- 7" Front rubber tyres
- Hard wearing nylon upholstery
- East to operate swing away leg rests
- Manual brake on rear wheels
- Robust design with 150kg user capacity
- Adjustable seat depth
- Available in Metallic Gun Metal Grey

## **Dimensions**

- Seat Width 18" (450mm) RG18TW
- Seat Width 20" (500mm) RG20TW
- Seat Depth 400mm 440mm
- Seat Height 500mm
  - Max User Capacity 150kg
- Wheelchair Weight 18" 19.5kg
- Wheelchair Weight 20" 21.5kg

## <u>Options Available</u>

- Anti tip wheels
- Stump support

Phone: 08 9248 4180 Fax: 08 9248 9452

E-mail: orders@redgumbrand.com.au

