

# ALPERTON

Dual motor lift lay-flat recliner with headrest and lumbar




Shown in Merino (Fresh)

Smart, compact and chic, the petite Alperton dual motor lift chair, with independent headrest and lumbar combines class with function. Available in two stylish colours, the Alperton takes comfort to a new level.

## Quick Information

- Independent headrest, legrest and lumbar function
- Footrest extension
- Steel seat box
- All close and lift button
- Convenient hand wand with USB connectivity

## Fabric Selection



**THEOREM**  
CONCEPTS

RECLINER AND LIFT CHAIR COLLECTION

## Specifications

Weight Capacity	158.7 kg
Seat Width	43.1 cm
Seat Depth	46.9 cm
Seat Height	44.4 cm
Chair Width	74.9 cm
Chair Depth	93.9 cm
Chair Height	44.4 cm
Fully Reclined Length	170 cm
Wall Space	40.6 cm
Foam Density	2.2
Cover	Merino
Cube	19.5
Box Size	76.9 cm x 76.9 cm x 80 cm
Net Weight	59 kg
Gross Weight	64.4 kg
Operation	USB Hand Wand